

# Bandage Series



## NINGBO HAOSHU HAORUN MEDICAL DRESSING CO.,LTD

ADD: NO1340 LIYUAN ROAD(N) HANSHU NINGBO ZHEJIANG CHINA 315000

Tel:0086-574-87706193

Fax:0086-574-87706193

Mobile:0086-13566542343

Email:sales@haorunmed.com

Website:www.haorunmed.com



NINGBO HAOSHU HAORUN MEDICAL DRESSING CO.,LTD



### Tubular bandage

Tubular bandage ( stockinette bandage ) is made of 100% cotton or with spandex , viscose other elastic fibers, soft and comfortable, different size is available. fit to feet. arm. elbow and so on.

#### Size width(natural) application

A	4.50cm	infant feet or arms
B	6.25cm	small hands and arms
C	6.75cm	medium arms small ankles
D	7.50cm	large arms .medium ankles. small kness
E	8.75cm	large ankles. medium knees. samll thighs
F	10.00cm	large knees, medium thighs
G	12.00cm	large thighs
J	17.50cm	small trunks
K	21.50cm	medium trunks
L	32.50cm	large trunks
M	37.50cm	extra-large trunks

We can supply specified size from clients.

Width : 1.5cm~33cm any length.

Packing: each roll in a OPP bag or box with print.



### Tubular bandage ( viscosa )

This kind of tubular bandage are soft .cool and elastic ,with a green line(or other color) in the middle of the bandage ,different size is available.

### Tubular bandage (latex free)

Tubular bandage latex free , with white color and skin color , different size is available .fit to feet. arm. elbow and so on.



### Tubular bandage (latex)

This kind of tubular bandage have a good elasticity ,plain and soft ,different size is available fit to feet, arm, elbow and so on.

### Stockinette

It is the best dressing product of tubular bandage ,high elasticity and good breathability .even in pressure ,tight to skin , different size is available.





**100% cotton elastic bandage**

100% cotton elastic bandage is the main product of our factory ,we are one of the largest manufacturers to produce the bandage in China

——materials: 100% high quality cotton

——color: natural White, bleached White, flesh, green, with or without blue/Red lines.

——weigh: 60–200g/m<sup>2</sup>

We can supply bandage with width ranging from 2.5cm to 30cm and length according to customers requirement .



**100% cotton plain elastic bandage ( raw edge )**

100% cotton plain elastic bandage have good elasticity ,good recovery ratio,good appearance and easily using. different size & color is available.



**100% cotton plain elastic bandage ( woven edge )**

100% cotton plain elastic bandage have good elasticity ,good recovery ratio,good appearance and easily using. different size & color is available.



**100% cotton crepe elastic bandage ( raw edge )**

100% cotton crepe elastic bandage have high absorbency and high elasticity, usually use in wound. surgical ,clinical, medical usage ,washable , different size & color is available.





#### Spandex elastic bandage

Spandex elastic bandage is made of natural cotton and spandex fiber with fixed ends, comfortable . for medical clinic .health care and athletic spots use .

- Width within 2.5–20cm and the length according to customs requirements
- Weight: 60–200g/m<sup>2</sup>
- Packing: each roll with clips, individual pack then a dozen in a bag or box, several bag/box in a carton
- Color: Nature White, bleached White flesh and so on



#### Spandex elastic bandage

spandex elastic bandage are made by cotton yarn and spandex threads , friendly to people ,to fix wound dressings for first aid. High elasticity ,washable ,sterilizable .



#### Spandex elastic bandage ( plain & crepe)

spandex elastic bandage are made by cotton yarn and spandex threads , friendly to people ,to fix wound dressings for first aid. High elasticity ,washable ,sterilizable.

### High elastic bandage

High elastic bandage is cool and comfortable to wear ,superior strength and elasticity ,with selvages and fixed ends.

- materials: cotton.polyester.rubber.
  - packing :each bandage is wrapped individuauy in waterproof bag,the outer package is strong cardboard carton to keep best storing condition .
  - weight : 60–200g/m<sup>2</sup>
  - different width and length available
- For more information, pls contract us



### High elastic bandage(latex)

Latex elastic bandage are orthopedic grade woven bandages that provide consistent elasticity,support and compression ,made of the finest quality yarns,these bandage can withstand numerous washing and sterillzation ,the material of latex rubber is gain authentication , that is innocuous to the skin of patient .



### High elastic bandage(latex free)

Latex free high elastic bandage ,orthopedic grade knitted bandage that can be wrapped and secured easily, it provides patient comfort and is durable enough to be reusable . For use in hospitals, physician office , school health and home care .



### PBT bandage

PBT bandage ( conforming gauze bandage ) , it applies to the orthopaedic surgeons, it feels comfortable 、lightness and elegant.good convenient, excellent elastic and reasonable price.

——It made of polyester fiber 、 polyamide 、 cotton or DTY

——Width within 2.5~20cm and the length according to customs requirements

packing: individual packing, then a dozen in a bag /box



### PBT elastic bandage ( thin )

PBT elastic bandage , plain or crepe . used to treat wound care or the daily accident .It's made of soft texture ,with high and permanent elasticity



### PBT elastic bandage ( thick )

PBT elastic bandage ( thick ) made of soft texture and cotton , with high and permanent elasticity ,suitable for the use of binding all parts of body, with woven edges, with smooth surface,200% elasticity ,it is ideal for fixation of wound dressing.



### Elastic bandage

This kind of bandage made of soft texture :spandex and cotton. Soft , elasticity 280~290% ,high absorbency, washable ,sterilizable,color can be white and flesh etc .

### Gauze Bandage

Gauze bandage (w o w bandage) cutting edges or woven edges , it has high absorbent and soft ,pure white ,widely used and cut out easily ,we can offer all kinds of gauze products with good quality and competitive price .

——material :pure absorbent cotton 21's .32's .40's

——mesh: 13 \17\20\24 threads



#### Gauze bandage ( woven edges )

side woven gauze bandage is made of 100% cotton with folded edges, soft and comfortable ,more convenient and more advantage .we can supply gauze bandage width within 1cm-20cm ,and any length according to customers .



#### Gauze bandage ( cut edges )

this gauze roll are made of pure cotton yarn ,through high temperature and pressure degreased and bleached ready-cut, it is the necessary products for hospital, White paper per roll .blue kraft per dozen .



### First aid bandage

First aid bandage usually used for emergency aid , it consists of absorbent pad and elastic bandage, soft and comfortable to the wound, easy to operate, various size for the bandage and the dressings for your choice.



### First aid bandage( thick)

First aid bandage unite the elastic bandage and the wound dressing, integrating the function of wrapping , stoping bleeding and accessorial fixing ,it can be used in different injures in different parts of body in case of battlefield or pre-hospital emergency.




### Plaster of paris bandage

Plaster of paris bandage is suitable for instant use , plaster covers the gauze completely, it is necessary product for surgery treatment.  
Width: 5cm/7.5cm/10cm/12.5cm/15cm/20cm, length: 2.7m/3.0m/4.6m



### Cohesive elastic bandage

Cohesive elastic bandage , excellent pressure bandage, without cutting off circulation, self adhesive to itself, no clips or fasteners needed, easily hand tearable. its made of nonwoven or cotton .different color & size is available.



### First aid bandage( thin)

Usually made of PBT bandage & absorbent pad. There is a thin membrane-like material in non-ticks pad on the outer side toward the wound ,it allows the body fluid to fully pass to the absorbent layer ,it also makes it easier to take the dressing from the wound . Available in multiple package, non sterile and sterile , individual pack in cellophane or pouch .



### Orthopaedic padding

Orthopaedic padding is wide-used supporting material for surgery and wadding of plaster of paris bandages, it is soft ,smooth and comfortable to body's shape .its made of cotton and polyester,  
Width: 5cm/7.5cm/10cm/12.5cm/15cm/20cm, length: 2.7m/3.0m/4.6m